



Exceeding Your  
Expectations  
We Build  
Trust

深圳市罗杰科技有限公司  
SHENZHEN SCIFLOW TECHNOLOGIES CO LTD  
[www.sciflow.com](http://www.sciflow.com) [eleflow@163.com](mailto:eleflow@163.com)

No Cleaning — No Cross Contamination — Reduce Costs

## Powder Spatulas

Easy to use spatulas - made in a cleanroom

New



- Made from white Polystyrene (FDA grade - conforms to CFR 177.1640)
- Tough, efficient 'U' shaped blade
- Individually wrapped
- Manufactured and packed in a medical cleanroom (Class 100,000) to ensure cleanliness
- Range of sizes
- Available pre-sterilised and non sterile
- Smooth ergonomic round handle



Part No.	Blade Length	Overall Length	Pre-Sterilised	Quantity per Box
8095A-150	150mm	230mm	No	100
8095A-250	250mm	330mm	No	100
8095A-150S	150mm	230mm	Yes	100
8095A-250S	250mm	330mm	Yes	100

## The Powder Lock System

The perfect combination for efficient surface sampling

New



- Powder Spatula complete with a quick connect Container
- Specially designed container to eliminate messy powder transfers from sampler to bottle
- Made from white Polystyrene (FDA grade - conforms to CFR 177.1640)
- Individually wrapped
- Manufactured and packed in a medical cleanroom (Class 100,000)
- Range of sizes
- Available pre-sterilised and non sterile
- Smooth ergonomic rounded handle
- Container is designed for easy, controlled powder discharge



Part No.	Blade Length	Overall Length	Pre-Sterilised	Quantity per Box
8097A-150	150mm	230mm	No	100
8097A-250	250mm	330mm	No	100
8097A-150S	150mm	230mm	Yes	100
8097A-250S	250mm	330mm	Yes	100


Easily controlled powder discharge



Powder Spatula & Container



Powder Spatula locked inside the container



Quick Clip for fast and secure connection - easy to undo when wearing gloves