



Exceeding Your  
Expectations  
We Build  
Trust

深圳市罗杰科技有限公司  
SHENZHEN SCIFLOW TECHNOLOGIES CO LTD  
[www.sciflow.com](http://www.sciflow.com) [eleflow@163.com](mailto:eleflow@163.com)

**No Cleaning — No Cross Contamination — Reduce Costs**

# Pipettes & Jumbo Pipettes




## Low Cost Disposable Liquid Samplers

- Two styles available - Thumb Tube Pipettes & Jumbo Pipettes
- Simple to use
- Individually wrapped
- Assembled in a medical cleanroom
- Available pre-sterilised

### Thumb Tubes or Jumbo Pipettes?

**Thumb Tubes:** Ideal for low volume sampling. This sampler is more likely to drip than the Jumbo Pipette

**Jumbo Pipettes:** The Jumbo Pipettes has tapered ends. This allows a large volume sample to be taken and ensures the operator has complete control. Also minimises/eliminates dripping.



### Thumb Tube Pipettes - small volume



Part No.	Overall Length	Overall Diameter	Internal Diameter	Pre-Sterilised	Material	Max Sample Volume*	Quantity per Box
5055A-50	1000mm	12mm	9mm	No	HDPE	60ml*	50
5055A-50S	1000mm	12mm	9mm	Yes	HDPE	60ml*	50

### Jumbo Pipettes - large volume, with extra control



Part No.	Overall Length	Overall Diameter	Internal Diameter	Pre-Sterilised	Material	Max. sample Volume*	Quantity per Box
8051H-450	450mm	21mm	17mm	No	HDPE	80ml*	50
8051H-450S	450mm	21mm	17mm	Yes	HDPE	80ml*	50
8051H-900	900mm	21mm	17mm	No	HDPE	180ml*	50
8051H-900S	900mm	21mm	17mm	Yes	HDPE	180ml*	50

\* The maximum sample volume will depend upon how deeply the sampler is inserted into the product